


Extend® Absorbable Synthetic Implants


Mid-Term Occlusion

Extend Absorbable Synthetic Implants provide occlusion for up to 3 months for patients experiencing short-term dry eye symptoms. These implants are ideal for the transient dry eye associated with LASIK or other ocular surgeries. Extend implants can also be used for allergy sufferers, contact lens discomfort, climate change or to sustain drug contact time and effectiveness.

A Simple Solution.

- Simple insertion technique
- No insertion tool needed; just forceps
- Can be used to trial permanent occlusion
- No foreign body sensation
- Made of a dissolvable copolymer of glycolic acid and trimethylene carbonate
- Billed/reimbursed the same as permanent silicone plugs



Extend Synthetic Absorbable Implants	Item #	Quantity
	0.2 mm x 2.0 mm	A13-202
	0.3 mm x 2.0 mm	A13-203
	0.4 mm x 2.0 mm	A13-204
	0.5 mm x 2.0 mm	A13-205
Absorbable Synthetic Implants (Variety Pack)	A13-200	10 pair box of: 3 pair — .3 mm 5 pair — .4 mm 2 pair — .5 mm

How To Order: Call 866-906-8080 or www.beaver-visitec.com

Beaver-Visitec International, Inc.
411 Waverley Oaks Road, Waltham, MA 02452
Tel: 1.866.906.8080 Fax: 1.866.906.4304
www.beaver-visitec.com

Beaver-Visitec International, Limited
85c Park Drive, Milton Park,
Abingdon, Oxfordshire, OX14 4RY, UK
Tel: 44.1865.601.256 Fax: 44.1865.595.761

 **Beaver Visitec**