

Plus™ Punctal Occluder

Designed to be of value to you and your patients




The Plus design is the result of BVI's proprietary micro-molding manufacturing process using a soft medical grade silicone. Fifteen years of clinical implantation reveal the design features of the Plus provide:

- A sharply angled nose for easy insertion
- Delicate contoured edges
- An ultra-smooth surface finish
- Excellent retention



Size Comparison Chart



| Plus Sterile Pre-Loaded Punctal Occluders | | Item # | Description |
|--|-----------|---------------------------|-------------|
|  <p>Occluders are packaged in a convenient (1) pair tray. Each occluder is individually sterile and pre-loaded on a high-quality disposable inserter.</p> | A14-204 | Small (0.4 mm - 0.55 mm) | |
| | A14-205 | Medium (0.55 mm - 0.7 mm) | |
| | A14-206 | Large (0.7 mm - 0.85 mm) | |
| | A14-207 | X-Large (0.85 mm +) | |
| Plus Bulk, Non-Sterile Punctal Occluders | | Item # | Description |
|  <p>Occluders are supplied non-sterile, 2 plugs per poly bag. 10 pair per box. 20 occluders total. Reusable inserters sold separately.</p> | A14-204BN | Small (0.4 mm - 0.55 mm) | |
| | A14-205BN | Medium (0.55 mm - 0.7 mm) | |
| | A14-206BN | Large (0.7 mm - 0.85 mm) | |
| | A14-207BN | X-Large (0.85 mm +) | |
| Reusable Insertion Instrument (for use with bulk occluders) | | Item # | Description |
|  <p>Non-Sterile Reusable Insertion Instrument for use with bulk occluders.</p> | A14-300 | 1 pair | |

How To Order: Call 866-906-8080 or www.beaver-visitec.com

Beaver-Visitec International, Inc.
 411 Waverley Oaks Road, Waltham, MA 02452
 Tel: 1.866.906.8080 Fax: 1.866.906.4304
www.beaver-visitec.com

Beaver-Visitec International, Limited
 85c Park Drive, Milton Park,
 Abingdon, Oxfordshire, OX14 4RY, UK
 Tel: 44.1865.601.256 Fax: 44.1865.595.761

